

The ibidi product family is comprised of a variety of μ-Slides and μ-Dishes, which have all been designed for high-end microscopic analysis of fixed or living cells.

The glass bottom versions of the μ-Slides and μ-Dishes are especially designed for TIRF and single molecule applications. The μ-Slide 2 well is an array of 2 square fields where cells can be cultivated and, subsequently, investigated with microscopical methods. It is intended for the optimization of experimental parameters like antibody dilution, seeding density, or the most effective drug concentration.

## Material

The glass bottom version of the μ-Slides are made of a standard μ-Slide but with a glass coverslip bottom. It is not possible to detach the bottom. The μ-Slides are not autoclavable since they are temperature stable only up to 80°C / 175°F.

### Optical Properties ibidi glass bottom

|                        |   |
|------------------------|---|
| Refractive index $n_D$ | 1.523   |
| Abbe number            | 55  |
| Thickness              | No. 1.5H (selected quality<br>170 μm, ± 5 μm) |
| Material               | Schott borosilicate glass,<br>D 263M          |

Protein coatings increase direct cell growth of adherent cells. Specific coatings on glass are possible following this protocol:

- Prepare your coating solution according to the manufacturer's specifications or reference. Prepare your μ-Slide. Adjust the concentration to a coating area of 7.5 cm<sup>2</sup> and 1.5 ml.
- Apply 1.5 ml into the growth area. Make sure that the entire bottom is covered with liquid easily tilting or shaking the μ-Slide. Put on the lid and leave at room temperature for at least 30 minutes.
- Aspirate the solution and wash with the recommended protein dilution buffer. Optionally, let dry at room temperature. Attention, some coating proteins might degenerate when drying!

## Geometry

The μ-Slide 2 well glass bottom provides a standard slide format according to ISO 8037/1.

### Geometry of μ-Slide 2 well glass bottom

|                                       |                     |
|---------------------------------------|---------------------|
| Number of wells                       | 2                   |
| Dimensions of wells (w × l × h) in mm | 21.2 × 23.3 × 9.3   |
| Growth area per well                  | 4.8 cm <sup>2</sup> |
| Coating area per well                 | 7.5 cm <sup>2</sup> |
| Recommended filling volume per well   | 1.5 ml              |
| Total height with lid                 | 10.8 mm             |
| Bottom matches coverslip              | No. 1.5             |

## Surface and coating

The μ-Slide 2 well glass bottom is manufactured with an uncoated glass coverslip. Washing steps (e.g. with PBS) before cell seeding can remove glass dust which is advantageous for direct cell growth on the surface.

## Solvents for Fixation and Staining

Cells can be observed live or fixed directly in the μ-Slide preferably on an inverted microscope. The slide material is compatible to acids, alkalis, PFA, and silicone oil. Alcohols may be used for short term incubation (e.g. cell fixation). Acetone is not compatible. Further specifications can be found at [www.ibidi.com](http://www.ibidi.com).

For optimal results in fluorescence microscopy and storage of stained probes ibidi provides a mounting medium (50001) optimized for μ-Dishes and μ-Slides.

## Seeding cells

- Trypsinize and count cells as usual. Dilute the cell suspension to the desired concentration. Depending on your cell type, application of a 5-11 × 10<sup>4</sup> cells/ml suspension should result in a confluent layer within 2-3 days.
- Apply 1.5 ml cell suspension into each well of the μ-Slide. Avoid shaking as this will result in inhomogeneous distribution of the cells.

- Cover the slide with the supplied lid. Incubate at 37°C and 5% CO<sub>2</sub> as usual.

Undemanding cells can be left in their seeding medium for up to three days and grow to confluence there. However, best results might be achieved when the medium is changed every 1–2 days. Carefully aspirate the old medium and replace it by 1.5 ml/well fresh medium.

**Tip:**

As you may know from the 96 well plates, a bent meniscus at the air–liquid interphase in small open wells will destroy the phase contrast effect of your microscope image. Use the Ph+ version to overcome this disturbing effect.

**Immersion Oil**

When using oil immersion objectives, use only the immersion oils specified in the table. The use of a non-recommended oil could lead to the damage of the plastic material and the objective.

| Company | Product          | Ordering Number  |
|---------|------------------|------------------|
| ibidi   | Immersion Oil    | (ibidi) 50101    |
| Zeiss   | Immersol 518 F   | (Zeiss) 444960   |
| Zeiss   | Immersol W 2010  | (Zeiss) 444969   |
| Leica   | Immersion liquid | (Leica) 11513859 |

μ-Slide 2 well Family


The μ-Slide 2 well family is available as open well and as a Ph+ version. See table below for choosing your μ-Slide 2 well.

μ-Slide 2 well



| Ordering Number | Treatment or Coating    | Characteristics                         |
|-----------------|-------------------------|---|
| 80286           | ibiTreat, sterile       | hydrophilic, tissue culture treated     |
| 80282           | Collagen IV, sterile    | protein coating                         |
| 80283           | Fibronectin, sterile*   | protein coating                         |
| 80284           | Poly-L-Lysine, sterile  | biopolymer coating                      |
| 80285           | Poly-D-Lysine, sterile* | biopolymer coating                      |
| 80281           | uncoated, sterile       | hydrophobic                             |
| 80287           | glass bottom            | glass coverslip No. 1.5H (170 μm ±5 μm) |

μ-Slide 2 well Ph+



| Ordering Number | Treatment or Coating    | Characteristics                         |
|-----------------|-------------------------|---|
| 80296           | ibiTreat, sterile       | hydrophilic, tissue culture treated     |
| 80292           | Collagen IV, sterile    | protein coating                         |
| 80293           | Fibronectin, sterile*   | protein coating                         |
| 80294           | Poly-L-Lysine, sterile  | biopolymer coating                      |
| 80295           | Poly-D-Lysine, sterile* | biopolymer coating                      |
| 80291           | uncoated, sterile       | hydrophobic                             |
| 80297           | glass bottom            | glass coverslip No. 1.5H (170 μm ±5 μm) |

\* available on request only

**For research use only!**

Further technical specifications can be found at [www.ibidi.com](http://www.ibidi.com). For questions and suggestions please contact us by e-mail [info@ibidi.de](mailto:info@ibidi.de) or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany.  
© ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.