



# Perifusion Tubing Set


## Model No. PERI-TUB



Three-stop  
Tubing



Assembly  
example



The Perifusion Tubing Set is specifically designed for the Biorep Perifusion System. Dead-volume is a critical parameter in perfusion experiments; therefore, tube lengths must be kept consistent between channels. At Biorep, each perifusion tubing set is manufactured to meet these strict length requirements.


Save time and money by using pre-assembled tubing sets for each perifusion experiment.

Sold in sets for 8 or 12 channel perifusion experiments (non-sterile).

- Temp range: -50 to 74°C (-58 to 165°F) Disposable
- Ultra-smooth inner bore
- Excellent biocompatibility
- Non-cytotoxic
- Non-pyrogenic
- Non-hemolytic
- Low protein binding

- **Specifications**

PERI-TUB	
<b>Material</b>	Tygon®
<b>OD</b>	0.091" (2.31 mm)
<b>Wall thickness</b>	0.036" (0.91 mm)
<b>ID</b>	0.015" (0.38 mm)



Tubing set up